




THD N-ANO BASIC



Use

THD N-Ano BASIC anoscope is a SINGLE USE medical device intended for the analysis of the anorectal canal.

General Description

THD N-Ano BASIC is a handy and light medical device, that presents in the handle a suitable housing useful for supporting the THD PEN-LIGHT (THD880070) or THD OPTI-FIBER (THD880012) lighthing systems.

Available in two dimensions (N-ANO e N-ANO MAXI) each of which can have truncated or bevelled extremity.

Characteristics

- In medical grade polycarbonate
- Complies with the biocompatibility requirements established by ISO 10993
- Atraumatic insertion profile
- Extremely compact
- Weight 40 gr
- Class I CE device
- Individually packaged
- Latex Free

Advantages:

- Small size, light weight and easy to manoeuvre
- Minimally invasive, atraumatic



Technical data sheet


REF. 800301

REF. 800302

REF. 800303

REF. 800304

EN



		N-ano	N-ano MAXI
d	Internal lumen	18.5 mm	22.5 mm
D	Outer diameter	21.7 mm	25.8 mm
L	Operational length	65 mm	72 mm
Lf	Bevel length	52 mm	47 mm
Lmax	Device length	91 mm	91 mm

THD N-Ano Basic					
TRUNCATED			BEVELLED		
Code	Q.ty	Trade Name	Description	N° RDM	EMDN
800301 - 50	50	THD N-ANO BASIC (TR)	N-Ano Basic Truncated Anoscope	2385382	G02060301
800302 - 50	50	THD N-ANO BASIC (BE)	N-Ano Basic Bevelled Anoscope	2385447	G02060301
800303 - 50	50	THD N-ANO BASIC MAXI (TR)	N-Ano Basic MAXI Truncated Anoscope	2385449	G02060301
800304 - 50	50	THD N-ANO BASIC MAXI (BE)	N-Ano Basic MAXI Bevelled Anoscope	2385668	G02060301